




FusionFit



Spine



EMPOWERING
HUMAN FITNESS



www.medifitdevices.com

Mono Axial Screw

Diameter
 Ø 4.5, 5.5, 6.5 mm

Length
 25 to 50 mm (5.0 mm Increment)

Catalogue No.
 TT-441



Mono Reduction Screw

Diameter
 Ø 4.5, 5.5, 6.5 mm

Length
 25 to 50 mm (5.0 mm Increment)

Catalogue No.
 TT-442



Poly Axial Screw

Diameter
 Ø 4.5, 5.5, 6.5 mm

Length
 25 to 50 mm (5.0 mm Increment)

Catalogue No.
 TT-443



Poly Reduction Screw

Diameter
 Ø 4.5, 5.5, 6.5 mm

Length
 25 to 50 mm (5.0 mm Increment)

Catalogue No.
 TT-444



Connecting Rod

Diameter
 Ø 5.5 mm

Length
 50 to 100 mm (10.0 mm Increment)
 125 to 200 mm (25.0 mm Increment)

Catalogue No.
 TT-446



Transverse Connector

Catalogue No.
 TT-448

Connecting Rod

Diameter
 Ø 3.0 mm

Length
 60 mm

Catalogue No.
 TT-445



Crosslink Connector

Length
 30 to 34mm, 34 to 40mm, 40 to 50mm, 50 to 70mm

Catalogue No.
 TT-447

Anterior Cervical Plate (H Shape)

Length
 21 to 29 mm (2.0 mm Increment)
 32 to 65 mm (3.0 mm Increment)
 70 mm, 75 mm

Catalogue No.
 TT-452



Cervical Screw (For Cervical Plate H Shape)



Diameter
 Ø 4.0 mm

Length
 10 to 26 mm (2.0 Increment)

Catalogue No.
 TT-446

TLIF Cage

Diameter
 Ø 8 to 12mm

Length
 25 mm, 30 mm

Catalogue No.
 TT-460




PLIF Cage

Diameter
 Ø 7 to 11 mm

Length
 20 mm, 25 mm

Catalogue No.
 TT-458



Kidney Cage

Degree
4°, 8°, 11°

Length
8 to 12 mm

Catalogue No.
TT-459



Bullet Cage

Diameter
Ø 8 to 12 mm


Length
25 mm, 30 mm

Catalogue No.
TT-467

Cervical Square Cage

Diameter
Ø 6 to 12 mm

Catalogue No.
TT-462



Cervical Dio Cage

Diameter
Ø 14 mm, 16 mm

Length
5 to 10 mm

Catalogue No.
TT-463

Cervical Round Cage

Diameter
 Ø 5 to 7 mm

Catalogue No.
 TT-464



Cervical Cage with Screw

Diameter
 Ø 12 mm, 14 mm

Length
 5 to 9 mm

Catalogue No.
 TT-465



Cervical Screw (For Cervical Cage)

Diameter
 Ø 4.0 mm

Length
 10 to 22 mm (2.0 mm Increment)

Catalogue No.
 TT-466



Expandable Cage

Diameter
 Ø 12 to 18 mm

Length
 20 to 25 mm, 25 to 35 mm, 35 to 50 mm, 50 to 75 mm

Catalogue No.
 TT-457



Eco Cage (Round Hole)

Diameter
 Ø 10 to 18 mm

Length
 15 to 60 mm (5.0 Increment)

Catalogue No.
 TT-456



Mesh Cage

Diameter
 Ø 10 to 18 mm

Length
 15 to 60 mm (5.0 Increment)

Catalogue No.
 TT-461

Staple Washer




Catalogue No.
 TT-449




Spine Instrument Set Box





Instrument	Quantity
Reduction Against Wrench	1
T Handle	1
Bolts Holding Rod with Double Head	2
Screwdriver	4
Rotation Device with Double Head	1
Bolts Holding Rod	1
Rod Holding Forceps (Curved)	1
Anti-rotation Socket	1
Hexagon Screwdriver 2.5mm	1
Single Side Broken off Device	1
Tap 4.5/5.5/6.5	3
Probe (Straight, Curved)	2
Hexagon Screwdriver 3.5/4.0	2
Rod Pressuring Device	1
Developing Pin Inserter (4 pins)	1
Double Joint Rotation Forceps	1
Compression Forceps	1
Distraction Forceps	1
Rod Bending Forceps	1
Long Mono & Poly Screw Break off Forceps	1
Gun Type Persuader	1
Spine Awl (Straight)	1
Ratchet Handle	1
Box	1




Cervical Plate Implants & Instruments Set





SIS 124	Instrument	Quantity
SIS 124-001	Cervical Pin Retractor	1
SIS 124-002	Bone Curette	1
SIS 124-003	Plate Holding Forcep	1
SIS 124-004	Drill & Tap Sleeve Adjustable	1
SIS 124-005	Plate Holding Pin	1
SIS 124-006	A.C.P. Plate Bender	1
SIS 124-007	Pin Holder	1
SIS 124-008	Q.C. T-Handle	1
SIS 124-009	Self Holding Screw Driver For 4.0mm Locking Screw	1
SIS 124-0010	Q.C.Drill	1
SIS 124-0011	Q.C.Bone Tap	1
SIS 124-0012	Lock Screw Driver	1
SIS 124-0013	Box	1




Eco & Expandable Cage Implants & Instrument Set




SIS 126	Instrument	Quantity
SIS 126-001	Cage Holder	1
SIS 126-002	Cage Pusher (Curved)	1
SIS 126-003	Cage Pusher (Straight)	1
SIS 126-004	Cage Expansion Bar (Straight)	1
SIS 126-005	Cage Expansion Bar (Curved)	1
SIS 126-006	Cage Lock Screw Holder	1
SIS 126-007	Box with Silicon Base	1




Cage Implants & Instrument Set



SIS 125	Instrument	Quantity
SIS 125-001	Chisel Reamer 7mm	1
SIS 125-002	Chisel Reamer 8mm	1
SIS 125-003	Chisel Reamer 9mm	1
SIS 125-004	Chisel Reamer 10mm	1
SIS 125-005	Chisel Reamer 11mm	1
SIS 125-006	Chisel Reamer 12mm	1
SIS 125-007	Sizer 7mm	2
SIS 125-008	Sizer 8mm	2
SIS 125-009	Sizer 9mm	2
SIS 125-010	Sizer 10mm	2
SIS 125-011	Sizer 11mm	2
SIS 125-012	Sizer 12mm	2
SIS 125-013	Graft Punch	1
SIS 125-014	PLIF & Kidney Cage Holder	1
SIS 125-015	TLIF Cage Holder	1
SIS 125-015.01	Bullet Cage Holder	1
SIS 125-016	Curette Dynamic (Straight)	1
SIS 125-017	Curette Dynamic (Left - Right)	1
SIS 125-018	Curette Square (Straight)	1
SIS 125-019	Curette Square (Left & Right)	1
SIS 125-020	Cage Pusher (Straight)	1
SIS 125-021	Cage Pusher (Curved)	1
SIS 125-022	Q.C.T- Handle	1
SIS 125-023	Bone Graft Holder	1
SIS 125-024	Box with Silicon Base	1





 medifit

Medifit Devices LLP

 +91-96015 17744

 info@medifitdevices.com

 www.medifitdevices.com

 C-1209, Titanium City Center, Near Sachin Tower,
Anandnagar Road, Ahmedabad, India - 380 015