




FusionFit



Spine



EMPOWERING
HUMAN FITNESS



www.medifitdevices.com

Mono Axial Screw

Diameter
 Ø 4.5, 5.5, 6.5 mm

Length
 25 to 50 mm (5.0 mm Increment)

Catalogue No.
 TT-441



Mono Reduction Screw

Diameter
 Ø 4.5, 5.5, 6.5 mm

Length
 25 to 50 mm (5.0 mm Increment)

Catalogue No.
 TT-442



Poly Axial Screw

Diameter
 Ø 4.5, 5.5, 6.5 mm

Length
 25 to 50 mm (5.0 mm Increment)

Catalogue No.
 TT-443



Poly Reduction Screw

Diameter
 Ø 4.5, 5.5, 6.5 mm

Length
 25 to 50 mm (5.0 mm Increment)

Catalogue No.
 TT-444



Connecting Rod

Diameter
Ø 5.5 mm

Length
50 to 100 mm (10.0 mm Increment)
125 to 200 mm (25.0 mm Increment)

Catalogue No.
TT-446



Transverse Connector

Catalogue No.
TT-448

Connecting Rod

Diameter
Ø 3.0 mm

Length
60 mm

Catalogue No.
TT-445



Crosslink Connector

Length
30 to 34mm, 34 to 40mm, 40 to 50mm, 50 to 70mm

Catalogue No.
TT-447

Anterior Cervical Plate (H Shape)

Length
 21 to 29 mm (2.0 mm Increment)
 32 to 65 mm (3.0 mm Increment)
 70 mm, 75 mm

Catalogue No.
 TT-452



Cervical Screw (For Cervical Plate H Shape)



Diameter
 Ø 4.0 mm

Length
 10 to 26 mm (2.0 Increment)

Catalogue No.
 TT-446

TLIF Cage

Diameter
 Ø 8 to 12mm

Length
 25 mm, 30 mm

Catalogue No.
 TT-460




PLIF Cage

Diameter
 Ø 7 to 11 mm

Length
 20 mm, 25 mm

Catalogue No.
 TT-458



Kidney Cage

Degree
4°, 8°, 11°

Length
8 to 12 mm

Catalogue No.
TT-459



Bullet Cage

Diameter
Ø 8 to 12 mm


Length
25 mm, 30 mm

Catalogue No.
TT-467

Cervical Square Cage

Diameter
Ø 6 to 12 mm

Catalogue No.
TT-462



Cervical Dio Cage

Diameter
Ø 14 mm, 16 mm

Length
5 to 10 mm

Catalogue No.
TT-463

Cervical Round Cage

Diameter
 Ø 5 to 7 mm

Catalogue No.
 TT-464



Cervical Cage with Screw

Diameter
 Ø 12 mm, 14 mm

Length
 5 to 9 mm

Catalogue No.
 TT-465

Cervical Screw (For Cervical Cage)

Diameter
 Ø 4.0 mm

Length
 10 to 22 mm (2.0 mm Increment)

Catalogue No.
 TT-466



Expandable Cage

Diameter
 Ø 12 to 18 mm

Length
 20 to 25 mm, 25 to 35 mm, 35 to 50 mm, 50 to 75 mm

Catalogue No.
 TT-457

Eco Cage (Round Hole)

- | Diameter
| Ø 10 to 18 mm
- | Length
| 15 to 60 mm (5.0 Increment)
- | Catalogue No.
| TT-456



Mesh Cage

- | Diameter
| Ø 10 to 18 mm
- | Length
| 15 to 60 mm (5.0 Increment)
- | Catalogue No.
| TT-461

Staple Washer




- | Catalogue No.
| TT-449




Spine Instrument Set Box





| Instrument | Quantity |
|--|----------|
| Reduction Against Wrench | 1 |
| T Handle | 1 |
| Bolts Holding Rod with Double Head | 2 |
| Screwdriver | 4 |
| Rotation Device with Double Head | 1 |
| Bolts Holding Rod | 1 |
| Rod Holding Forceps (Curved) | 1 |
| Anti-rotation Socket | 1 |
| Hexagon Screwdriver 2.5mm | 1 |
| Single Side Broken off Device | 1 |
| Tap 4.5/5.5/6.5 | 3 |
| Probe (Straight, Curved) | 2 |
| Hexagon Screwdriver 3.5/4.0 | 2 |
| Rod Pressuring Device | 1 |
| Developing Pin Inserter (4 pins) | 1 |
| Double Joint Rotation Forceps | 1 |
| Compression Forceps | 1 |
| Distraction Forceps | 1 |
| Rod Bending Forceps | 1 |
| Long Mono & Poly Screw Break off Forceps | 1 |
| Gun Type Persuader | 1 |
| Spine Awl (Straight) | 1 |
| Ratchet Handle | 1 |
| Box | 1 |




Cervical Plate Implants & Instruments Set





| SIS 124 | Instrument | Quantity |
|--------------|---|----------|
| SIS 124-001 | Cervical Pin Retractor | 1 |
| SIS 124-002 | Bone Curette | 1 |
| SIS 124-003 | Plate Holding Forcep | 1 |
| SIS 124-004 | Drill & Tap Sleeve Adjustable | 1 |
| SIS 124-005 | Plate Holding Pin | 1 |
| SIS 124-006 | A.C.P. Plate Bender | 1 |
| SIS 124-007 | Pin Holder | 1 |
| SIS 124-008 | Q.C. T-Handle | 1 |
| SIS 124-009 | Self Holding Screw Driver For 4.0mm Locking Screw | 1 |
| SIS 124-0010 | Q.C.Drill | 1 |
| SIS 124-0011 | Q.C.Bone Tap | 1 |
| SIS 124-0012 | Lock Screw Driver | 1 |
| SIS 124-0013 | Box | 1 |




Eco & Expandable Cage Implants & Instrument Set




| SIS 126 | Instrument | Quantity |
|-------------|-------------------------------|----------|
| SIS 126-001 | Cage Holder | 1 |
| SIS 126-002 | Cage Pusher (Curved) | 1 |
| SIS 126-003 | Cage Pusher (Straight) | 1 |
| SIS 126-004 | Cage Expansion Bar (Straight) | 1 |
| SIS 126-005 | Cage Expansion Bar (Curved) | 1 |
| SIS 126-006 | Cage Lock Screw Holder | 1 |
| SIS 126-007 | Box with Silicon Base | 1 |




Cage Implants & Instrument Set



| SIS 125 | Instrument | Quantity |
|----------------|--------------------------------|----------|
| SIS 125-001 | Chisel Reamer 7mm | 1 |
| SIS 125-002 | Chisel Reamer 8mm | 1 |
| SIS 125-003 | Chisel Reamer 9mm | 1 |
| SIS 125-004 | Chisel Reamer 10mm | 1 |
| SIS 125-005 | Chisel Reamer 11mm | 1 |
| SIS 125-006 | Chisel Reamer 12mm | 1 |
| SIS 125-007 | Sizer 7mm | 2 |
| SIS 125-008 | Sizer 8mm | 2 |
| SIS 125-009 | Sizer 9mm | 2 |
| SIS 125-010 | Sizer 10mm | 2 |
| SIS 125-011 | Sizer 11mm | 2 |
| SIS 125-012 | Sizer 12mm | 2 |
| SIS 125-013 | Graft Punch | 1 |
| SIS 125-014 | PLIF & Kidney Cage Holder | 1 |
| SIS 125-015 | TLIF Cage Holder | 1 |
| SIS 125-015.01 | Bullet Cage Holder | 1 |
| SIS 125-016 | Curette Dynamic (Straight) | 1 |
| SIS 125-017 | Curette Dynamic (Left - Right) | 1 |
| SIS 125-018 | Curette Square (Straight) | 1 |
| SIS 125-019 | Curette Square (Left & Right) | 1 |
| SIS 125-020 | Cage Pusher (Straight) | 1 |
| SIS 125-021 | Cage Pusher (Curved) | 1 |
| SIS 125-022 | Q.C.T- Handle | 1 |
| SIS 125-023 | Bone Graft Holder | 1 |
| SIS 125-024 | Box with Silicon Base | 1 |






 medifit

Medifit Devices LLP

 +91-96015 17744

 info@medifitdevices.com

 www.medifitdevices.com

 C-1209, Titanium City Center, Near Sachin Tower,
Anandnagar Road, Ahmedabad, India - 380 015